

Revised Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® PGT-A Trisomy 18 Reference Material

Material Number: 0720-0776

Batch Number: 10683803

Material Description: Genomic DNA extracted from a male human cell line carrying three copies of chromosome 18.

Concentration (Qubit dsDNA HS Assay): 23.4 pg/μL

Fill Volume: 10 μL

Date of Manufacture: 23 OCT 2023


Expiration Date: 24 AUG 2026

Trisomy Call using commercially available NGS-based PGT-A: Trisomy 18

Sex determined using commercially available NGS-based PGT-A: Male

Storage: -20°C

Approval:

A handwritten signature in black ink, appearing to be "M H".

24 SEP 2024

Prepared By

Date