

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq[®] PGT-A Trisomy 18 Reference Material

Material Number: 0720-0776 Batch Number: 10613698

Material Description: Genomic DNA extracted from a male human cell line carrying three copies of chromosome 18.

Concentration (Qubit dsDNA HS Assay): 32.2 pg/μL

Fill Volume: 10 μL


Date of Manufacture: 26 APR 2022 Expiration Date: 24 AUG 2023

Trisomy Call using commercially available NGS-based PGT-A: Trisomy 18

Sex determined using commercially available NGS-based PGT-A: Male

Storage: -20°C

Approval:



Prepared By

26 APR 2022

Date