

Revised Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® PGT-A Trisomy 13 Reference Material

Material Number: 0720-0777

Batch Number: 10683804

Material Description: Genomic DNA extracted from a male human cell line carrying three copies of chromosome 13.

Concentration (Qubit dsDNA HS Assay): 27.4 pg/μL

Fill Volume: 10 μL

Date of Manufacture: 23 OCT 2023

Expiration Date: 25 AUG 2026

Trisomy Call using commercially available NGS-based PGT-A: Trisomy 13

Sex determined using commercially available NGS-based PGT-A: Male

Storage: -20 °C

Approval:

A handwritten signature in black ink, appearing to be "M H".

24 SEP 2024

Prepared By

Date