

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® PGT-A Euploid Reference Material

Material Number: 0720-0778 Batch Number: 10613700

Material Description: Genomic DNA extracted from a normal female human cell line.

Concentration (Qubit dsDNA HS Assay): 33.6 pg/μL

Fill Volume: 10 μL

Date of Manufacture: 27 APR 2022 Expiration Date: 25 AUG 2023

Trisomy Call using commercially available NGS-based PGT-A: Euploid

Sex determined using commercially available NGS-based PGT-A: Female

Storage: -20°C

Approval:

A handwritten signature in blue ink, appearing to be "M. J.", written over a horizontal line.

28 APR 2022

Prepared By

Date