

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 13 Male - Matched Reference Material
Material Number: 0720-0779 Batch Number: 10551201
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 34.3 ng/mL
Volume: 1.0 mL
Date of Manufacture: 19 APR 2021 Expiration Date: 13 MAY 2024
Fetal Fraction by NGS (Illumina VeriSeq): 12.2%

Trisomy Call for Chromosome 13 by NGS (Illumina VeriSeq): Pass
Call for fetal gender by NGS (Illumina VeriSeq): XY
Storage: 2-8 °C

Approval:

A handwritten signature in blue ink, appearing to be "J. M. M." or similar.

26 APR 2021

Prepared By

Date