

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Trisomy 13 Male - Matched Reference Material  
Material Number: 0720-0779                      Batch Number: 10766380  
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).  
Concentration (Qubit dsDNA BR Assay): 33.2 ng/mL  
Volume: 1.0 mL  
Date of Manufacture: 04 JUN 2025                      Expiration Date: 16 APR 2028  
Fetal Fraction by NGS (Chroma Panel) 13%  
  
Trisomy Call for Chromosome 13 by NGS (Chroma Panel) Pass  
Call for fetal gender by NGS (Chroma Panel) XY  
Storage: 2-8 °C

**Approval:**

A handwritten signature in blue ink, appearing to be "DP".

02 JUL 2025

Prepared By

Date