

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 13 Male - Matched Reference Material
Material Number: 0720-0779 Batch Number: 10788501
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 31.1 ng/mL
Volume: 1.0 mL
Date of Manufacture: 16 OCT 2025 Expiration Date: 16 APR 2028
Fetal Fraction by NGS (Chroma Panel): 13%
Trisomy Call for Chromosome 13 by NGS (Chroma Panel): Pass
Call for fetal gender by NGS (Chroma Panel): XY
Storage: 2-8 °C

Approval:

A handwritten signature in blue ink, appearing to read "Brianne Dwyer".

Prepared By

27 Oct 2025

Date