

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 13 Male - Matched Reference Material
Material Number: 0720-0779 Batch Number: 10810848
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 27.1 ng/mL
Volume: 1.0 mL
Date of Manufacture: 24 APR 2026 Expiration Date: 05 MAR 2030
Fetal Fraction by NGS (Chroma Panel): 16%
Trisomy Call for Chromosome 13 by NGS (Chroma Panel): Pass
Call for fetal gender by NGS (Chroma Panel): XY
Storage: 2-8 °C

Approval:

A handwritten signature in blue ink, appearing to be "M. H.", written over a horizontal line.

29 APR 2026

Prepared By

Date