

# Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Triple X Syndrome Reference Material

Material Number: 0720-1069 Batch Number: 10681848

Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma)

Date of Manufacture: 15 AUG 2023 Expiration Date: 25 JUL 2027

Fill Volume: 1.0 mL

Storage: 2-8 °C

Concentration (Invitrogen™ Qubit™ dsDNA BR Assay): 31.9 ng/mL

Method Used for NIPT: Illumina VeriSeq™ NIPT Solution v2

Fetal Fraction: 21.0 %

Call for Triple X: Pass

**Approval:**



17 AUG 2023

Prepared By

Date