

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Triple X Syndrome Reference Material

Material Number: 0720-1069

Batch Number: 10816831

Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma)

Date of Manufacture: 01 JUL 2026

Expiration Date: 25 JUL 2027

Fill Volume: 1.0 mL

Storage: 2-8 °C


Concentration
(Invitrogen™ Qubit™
dsDNA BR Assay): 31.9 ng/mL

Method Used for
NIPT: Illumina Veriseq™ NIPT Solution v2

Fetal Fraction: 21.0 %

Call for Triple X: Pass

Approval:

A handwritten signature in black ink, appearing to be "MJA", written over a horizontal line.

14 JUL 2026

Prepared By

Date