

Certificate of Analysis

Product Description: HAV IgM Plasma
 Material No(s): 0301-0006, 0301-0007, 0301-0008
 Lot No: 10128180
 Donor ID: NA
 Gender: NA
 Age: NA
 Country of Origin: Germany
 Collection Date: 21-Apr-10
 Anticoagulant: 4% Sodium Citrate
 Storage Temperature: -20 degrees C
 Shipping Conditions: Dry ice
 Intended Use: For research use only. Not for use in diagnostic procedures.
 Safety Precautions: Use universal precautions when handling this material.
 Test Method(s): See below

Lot Analysis					
Test Date	Analyte	Test Method/Manufacturer	Result	Units	Notes
12-Aug-15	HAV IgM	DIASORIN EIA	5.186	S/CO	NA
This product was tested for anti-HIV 1/2, anti-HCV, and HBsAg and found to be negative.					

Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.
 NA = Not Available

Approval:

Jessica Calae

10-Sep-15

Prepared By

Date

Christine Pawia

11-Sep-15

QA Verified By

Date

