

Certificate of Analysis

Product Description: HIV-1 Subtype C
 Material No(s): 0315-0051, 0315-0052, 0315-0053
 Lot No: 9145671
 Donor ID: NA
 Demographics: NA
 Collection Date: 19-Aug-2002
 Anticoagulant: NA
 Storage Temperature: -70 degrees C
 Shipping Conditions: Dry ice
 Intended Use: For research use only. Not for use in diagnostic procedures.
 Safety Precautions: Use universal precautions when handling this material.
 Test Method(s): See below

| Lot Analysis | | | | |
|---------------|------------|--|--------|-------|
| Test Date | Analyte | Test Method/Manufacturer | Result | Units |
| 23-April-2010 | HIV-1 GENO | Reference Lab (Siemens), Tru-Geno Method | C | NA |
| 19-Oct-2010 | HIV-1 RNA | Roche COBAS AmpliPrep TaqMan v1.0 | 18700 | C/mL |
| 15-Nov-2010 | HIVC | Bio-Rad EIA a-HIV 1/2 | 10.381 | S/CO |

This product was tested for HBSAG, and HCV and found to be negative.

Data are provided for informational purposes.
 SeraCare Life Sciences does not claim that others can duplicate test results exactly.
 NA = Not Available

Approval:




QA _____ Date

