

## Certificate of Analysis

Product Description:	HIV-1 Subtype A Plasma
Material No(s):	0315-0045, 0315-0046, 0315-0047
Lot No:	9250622
Donor ID:	BD110650
Demographics:	Kyrgyzstan
Collection Date:	09-FEB-11
Anticoagulant:	Heparin
Storage Temperature:	-70 degrees C
Shipping Conditions:	Dry ice
Intended Use:	For research use only. Not for use in diagnostic procedures.
Safety Precautions:	Use universal precautions when handling this material.
Test Method(s):	See below

Lot Analysis				
Test Date	Analyte	Test Method/Manufacturer	Result	Units
16-Mar-11	HIV-1 RNA	Roche COBAS AmpliPrep TaqMan 1.0	50900	C/mL
21-Mar-11	a-HIV 1/2	Bio-Rad EIA	11.07	S/CO

This product was tested for HBsAg, and a-HCV and found to be negative.

Data are provided for informational purposes.  
 SeraCare Life Sciences does not claim that others can duplicate test results exactly.  
 NA = Not Available

**Approval:**

 _____ QA	24 NOV 11 _____ Date
--	----------------------------

