

Certificate of Analysis

Revision 02

Product Description: dsDNA Antibody Positive Plasma
 Material No(s): 0340-0029, 0340-0030, 0340-0031, 0340-0032
 Lot No: 9260435
 Donor ID: BD111203
 Gender: NA
 Age: NA
 Country of Origin: NA
 Collection Date: 29-Jun-12
 Anticoagulant: 4% Sodium Citrate
 Storage Temperature: -20 degrees C
 Shipping Conditions: Dry ice
 Intended Use: For research use only. Not for use in diagnostic procedures.
 Safety Precautions: Use universal precautions when handling this material.
 Test Method(s): See below

Lot Analysis					
Test Date	Analyte	Test Method/Manufacturer	Result	Units	Notes
06-Aug-12	dsDNA Antibody	Diamedix EIA	> 200	IU/mL	ALR
06-Aug-12	dsDNA Antibody	Diamedix EIA	82.24	IU/mL	1:10 dilution
This product was tested for HBsAg, Antibody to HIV 1/2, Antibody to HCV, HIV1 RNA, HBV DNA, and HCV RNA, and found to be negative. The donor has reported a diagnosis of SLE .					

Additional Data					
Test Date	Analyte	Test Method/Manufacturer	Result	Units	Notes
04-Feb-13	C Reactive Protein	Beckman Coulter Immage	0.566	mg/dL	

Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.
 NA = Not Available

Approval:

Lynne Labonte _____ 14 Dec 15
 Prepared By _____ Date
Tsering Chelung _____ 14 Dec 15
 QA Verified By _____ Date

